

Segmenting Large Electron Microscopic Image Volumes

An Introduction to NBCR image analysis and segmentation tools



National Biomedical Computation Resource Summer Training Program @ UC San Diego

Chris Churas & Matthias Haberl (original instructions created by Alex Perez P.H.D)

August 7th, 2017

churas@ncmir.ucsd.edu & mhaberl@ucsd.edu

Hands-On Session 1

The Manual and Automated Segmentation of Organelles in 3D EM Data


Goals

By the end of this session, you will be able to:

1. Use IMOD program to get information about SBEM image stack
2. Generate a training dataset from training data with CHMutil/IMOD that can be run by CHM


1. Download and Unzip the Dataset

- A. Download the ZIP file entitled *download_dataset_1.zip* from the BBDTC portal
- B. Unzip the contents to an easily accessible location:
 - i. For Windows 7/8: Unzip the file into a new folder in your home directory. If you used the IMOD cygwin install, the home directory will be C:\Users\<username> (for windows 10 it may be C:\cygwin\home\<username>)
 - ii. For Mac: Unzip the file into a new folder (e.g. /Users/<username>/download_dataset_1).
 - iii. For Linux: Unzip the file into a new folder (e.g. /home/<username>/download_dataset_1).
- C. Open your terminal program
- D. Enter the directory to which you unzipped the dataset by using the *cd* command in the terminal:
cd ~/download_dataset_1
- E. List the directory's contents using the *ls* command. You should see one file, *sbem_train.mrc*



```
smithy@smith ~  
$ ls Downloads/download_dataset_1.zip  
Downloads/download_dataset_1.zip  
  
smithy@smith ~  
$ unzip Downloads/download_dataset_1.zip  
Archive: Downloads/download_dataset_1.zip  
  creating: download_dataset_1/  
    inflating: download_dataset_1/sbem_train.mrc  
  
smithy@smith ~  
$ cd download_dataset_1/  
  
smithy@smith ~/download_dataset_1  
$ ls  
sbem_train.mrc  
  
smithy@smith ~/download_dataset_1  
$ |
```

Example Output
Windows 7/8



```
smith@DESKTOP-ECJR8SM ~  
$ unzip $HOMEPATH/Downloads/download_dataset_1.zip  
Archive: \\Users\smith\Downloads\download_dataset_1.zip  
  creating: download_dataset_1/  
    inflating: download_dataset_1/sbem_train.mrc  
  
smith@DESKTOP-ECJR8SM ~  
$ cd download_dataset_1/  
  
smith@DESKTOP-ECJR8SM ~/download_dataset_1  
$ ls sbem_train.mrc  
sbem_train.mrc  
  
smith@DESKTOP-ECJR8SM ~/download_dataset_1  
$ |
```

Example Output
Windows 10

These screen shots show use of *unzip* command from terminal to unzip the dataset file on Windows 7/8/10

2. Determine the Properties of the Dataset

- A. At the command line, type `header sbem_train.mrc`
- The program `header` is part of the IMOD software suite*, and reads pertinent details from MRC files. The output should contain the following:

```

Number of columns, rows, sections ..... 500 500 2
Map mode ..... 0 (bytes - signed in file)
Start cols, rows, sects, grid x,y,z ... 0 0 0 500 500
2
Pixel spacing (Angstroms)..... 77.98 77.98 300.0
Cell angles ..... 90.000 90.000 90.000
Fast, medium, slow axes ..... X Y Z
Origin on x,y,z ..... 0.000 0.000 0.000
Minimum density ..... 119.00 ( -11.000 in file)
Maximum density ..... 180.00 ( 54.000 in file)
Mean density ..... 152.41 ( 24.446 in file)
tilt angles (original,current) ..... 0.0 0.0 0.0 0.0 0.0 0.0
Space group,# extra bytes,idtype,lens . 0 0 0 0

```


Voxel dimensions of the file, in (X Y Z).

0 = 8-bit, unsigned integers
1 = 16-bit, signed integers
2 = 32-bit, floating point
6 = 16-bit, unsigned integers

Voxel sizes, in Angstroms per voxel, in (X Y Z). This dataset has anisotropic voxels, with 7.798 nm/voxel in X and Y, and 30 nm/voxel in Z. Anisotropic voxels are common in SBEM datasets.

* For an exhaustive list of IMOD programs and help documents, visit: http://bio3d.colorado.edu/imod/doc/program_listing.html

Directory Structure:



```

~/download_dataset_1
Example Output
smithy@smith ~/download_dataset_1
$ header sbem_train.mrc

RO image file on unit 1 : sbem_train.mrc      Size=      489 K

Number of columns, rows, sections .....      500      500      2
Map mode .....      0 (byte)
Start cols, rows, sects, grid x,y,z ...      0      0      0      500      500      2
Pixel spacing (Angstroms).....      77.98      77.98      300.0
Cell angles .....      90.000      90.000      90.000
Fast, medium, slow axes .....      X      Y      Z
Origin on x,y,z .....      -999.0      0.000      -14.00
Minimum density .....      119.00
Maximum density .....      180.00
Mean density .....      152.41
tilt angles (original,current) .....      0.0 0.0 0.0 0.0 0.0 0.0
Space group,# extra bytes,idtype,lens .      0      0      0      0

10 Titles :
raw2mrc: Converted to mrc format.                  19-Feb-12 12:22:31
NEWSTACK: Images copied                             21-Feb-12 11:37:20
NEWSTACK: Images copied                             21-Feb-12 12:19:51
NEWSTACK: Images copied                             22-Feb-12 08:54:59
NEWSTACK: Images copied                             24-Feb-12 12:04:59
NEWSTACK: Images copied, transformed, floated to means 28-Feb-12 22:16:28
BOXSTARTEND: 10 starts clipped out and averaged      3-Aug-15 13:50:47
NEWSTACK: Images copied                             3-Aug-15 13:52:31
NEWSTACK: Images copied                             3-Aug-15 13:52:50
NEWSTACK: Images copied                             14-Jul-17 08:47:41

smithy@smith ~/download_dataset_1
$

```

3. Create IMOD Model File for Training Labels

- A. Open the training data stack in 3dmod: *3dmod sbem_train.mrc*
- B. In this course, we will be generating training data for mitochondria. Using the techniques described in the preceding lecture, manually segment all instances of mitochondria in the two tiles of your training data stack.

As an additional reference, you can view the file:

Generating_CHM_training_data_with_IMOD.pdf



Or visit the following URLs:

<http://bio3d.colorado.edu/imod/doc/3dmodguide.html> – Detailed information about manual segmentation and model file structure in IMOD.

<https://www.youtube.com/watch?v=BsNSVLIQ-cE> – Useful video illustrating the use of IMOD's Drawing Tools

- C. Save the model file as *sbem_train.mod*

Directory Structure:



4. Generate a Training Label Stack

The goal of this step is to generate a new MRC stack with the same dimensions as the training images we previously segmented. However, in this stack, all pixels inside of the traced contours will have values of one, and all pixels outside of the traced contours will have values of zero. In this way, a binary label stack will be created that will serve to tell the CHM training algorithm where mitochondria are. We will use the IMOD program *imodmop* to generate this stack.

- A. Create the binary label stack using the following command:

```
imodmop -mask 1 sbem_train.mod sbem_train.mrc sbem_train_labels.mrc
```



To understand the meaning of these arguments and the *imodmop* syntax, view the *imodmop* man page in the terminal by typing *man imodmop*

- B. Visualize the label stack using 3dmod:

```
3dmod sbem_train_labels.mrc
```

- C. In the 3dmod window, select Image → Pixel View. Right-click on some of the white blobs in the Zap window, and confirm that their pixels are one-valued.

Directory Structure:



5. Generate Training Image and Label PNGs

The CHM algorithm requires inputs in the form of PNG files. Thus, before we can continue to job submission with the SLASH portal, we must first convert from MRC to PNG.

- A. First, we need to make sub-directories under the *train* folder for images and labels:

```
mkdir -p train/images  
mkdir train/labels
```


- B. Next, convert the MRC stack of raw training images to individually numbered PNGs using the IMOD program *mrc2tif*. The '-p' argument forces conversion to PNG, rather than TIF and the 'x' at end is needed by *mrc2tif* as a filename prefix:

```
mrc2tif -p sbem_train.mrc train/images/x
```

- C. Convert the MRC stack of training labels to individually numbered PNGs using the IMOD program *mrc2tif*:

```
mrc2tif -p sbem_train_labels.mrc train/labels/x
```

Directory Structure:



```
~/download_dataset_1  
smithy@smith ~/download_dataset_1  
$ mkdir -p train/images  
smithy@smith ~/download_dataset_1  
$ mkdir train/labels  
smithy@smith ~/download_dataset_1  
$ mrc2tif -p sbem_train.mrc train/images/x  
Writing PNG images. ..  
smithy@smith ~/download_dataset_1  
$ mrc2tif -p sbem_train_labels.mrc train/labels/x  
Writing PNG images. ..  
smithy@smith ~/download_dataset_1  
$ ls train/images/  
x.000.png  x.001.png  
smithy@smith ~/download_dataset_1  
$ ls train/labels/  
x.000.png  x.001.png  
smithy@smith ~/download_dataset_1  
$ |
```